



Mobile Point of Care
10 inch ruggedized
mobile clinical assistant

CliniScope

FIMI



MCA104

A life time of care at the touch of a screen

The CliniScope mobile clinical assistant enables nurses and physicians to take more informed clinical decisions at the point of care by using an ultra-ergonomic and medical grade always connected touch-screen device that integrates all the features needed to increase patient safety and reduce administrative workload.

Ultra-ergonomic

- lightweight design (approx. 1,8 kg) with back handgrip to minimize fatigue and grab & go docking station
- 10 inch XGA pen-touch-screen combined with intuitive controls and LED indicators
- smart swappable battery for seamless continuous operation

All-in-one functionality

- Microsoft® Genuine Windows® XP Tablet PC Edition or Genuine Windows® Vista® Business operating system
- Intel® Core™ Solo processor (1,2GHz), 60GB HDD, 1 GB DDR2 RAM
- built-in barcode scanner and RFID reader to support the five rights of patient safety
- digital camera (2Mpixel) to capture pictures of patient's wounds
- Bluetooth® for cordless connection of external devices and vital sign sensors
- USB port behind rubber cover to interface with USB devices when not docked
- optional docking station with built-in battery charger, 3 USB ports and Ethernet connectivity

Medical appliance

- easy to wipe clean to reduce spread of infection
- medical grade certification for maximum patient protection
- durable design for less risk of downtime (1m drop proof)

PHILIPS
sense and simplicity

